

Kinase Assay

UKC™ 1.0 – Universal Kinase Continuous Assay
Real-time, standardized detection of kinase activity

Key Features

Universal Kinase Compatibility

Supports diverse kinase classes—serine/threonine, tyrosine, lipid, and metabolic kinases—independent of substrate identity.

Continuous Real-Time Readout

NADPH fluorescence is monitored continuously, enabling precise measurement of phosphorylation rates, catalytic efficiency, and substrate preference.

Robust Coupled-Enzyme System

An optimized combination of sugar kinases and glucose-6-phosphate dehydrogenase converts **ADP release** into a strong, highly sensitive fluorescent signal.

Standardized Platform for Inhibitor Discovery

Fixed reaction chemistry ensures high reproducibility, making UKC™ 1.0 suitable for **high-throughput inhibitor screening**, **IC₅₀ determination**, **mode-of-inhibition studies**, and comparative profiling of inhibitor or substrate libraries.

Assay Principle

The UKC™ assay couples the kinase phosphorylation reaction to NADPH formation.

As a kinase transfers a phosphate group from ATP to its substrate, **ADP is released**.

In the coupled system, ADP is rapidly converted through sugar kinase activity into intermediates that drive **NADP⁺ → NADPH** reduction.

The **rate of NADPH fluorescence increase directly reflects kinase activity**, enabling real-time kinetic analysis.

Kit Components

Enzyme Cocktail

- Sugar kinases optimized to efficiently convert ADP released from kinase reactions
- Glucose-6-phosphate dehydrogenase tuned for continuous NADP⁺ reduction and stable fluorescent output

Assay Buffer

- Formulated for broad kinase compatibility, including required intermediate substrates

Optional Add-Ons (upon request)

- Positive control kinase
- Positive control peptide or protein substrate

User Guide

- Detailed protocol for kinetic assays, substrate profiling, and high-throughput inhibitor screening (HTS)

